



QUANTUM®

The Rehab Power Chair Solutions Company

Q4

Exceptional features,
combined with a
compact chassis



Images illustrated:

1. Rear of the Q4
2. ATX Suspension
3. Manual elevating and swing-away legrest
4. Power Recline
5. Q4 with Power Tilt and power Recline with a 18"x18" ION seat

* images used for illustration purposes only



CRASH TESTED



The Rehab Power Chair Solutions Company

Tel: 01869 324600
Email: sales@quantumrehab.co.uk
Web: www.quantumrehab.co.uk



Tel: 01869 324600
Email: sales@quantumrehab.co.uk
Web: www.quantumrehab.co.uk

Introduction

The Quantum® Q4 is a compact 4mph Mid-Wheel Drive power chair, designed for indoor and outdoor use. The Q4 has been specifically designed with the new ION Seat, which features two individual seat sizes.

The Q4 features mid-wheel 6 design, which enables tight turning radius and outstanding performance making the Q4 an excellent choice for the demands of the most active client.

Seating Options

Featuring a soft yet flexible cushioning system the Quantum® Q4 is available in two seat sizes.

16" x 16"
40.64 cm x 40.64 cm

18" x 18"
45.72 cm x 45.72 cm

Both sizes feature their own unique seating systems (See Power and Manual Options).



Drive System and Suspension



Featuring a 2-motor, intuitive drive system along with the unique Mid-Wheel 6® drive system, which allows six wheels on the ground for maximum stability.

The Quantum® Q4 also comes as standard with the exceptional ATX Suspension system allowing for a smoother journey.

Power and Manual Options

Offering two separate Powered and Manual seating options (See below). Each configuration features its' own standard features unique to the user.



Image illustrates Manual Recline



Option 1:

Static Tilt | Manual Recline:

This option features an easy attendant recline control system, and comes as standard with swing-away legrests.



Option 2:

Power Tilt | Power Recline:

This option features a quick and easy to use powered positioning system, along with manual elevating legrests and headrest as standard.

Technical Specifications

Drive Wheels (cm) 35.56 cm pneumatic or solid	OMNI Casters (cm) Front and Rear: 15.25 cm	Max Speed (cm) 4 Mph (6.4 km/h)	Ground Clearance (cm) Centre of frame: 7.62 cm	Turning Radius (cm) 50.8 cm	Overall Length (cm) 90.17 cm	Overall Width (cm) 61 cm	Ion Seating Options (inches) 16" x 16" or 18" x 18"	Batteries x 2 volt 50 - 55 AH Deep Cycle	Available Electronics 70A PG VR2 4 or 6 key	Weight Capacity (stone) 21 Stone 6 lbs (137kg)
Base Weight (kg) 62 kg	Battery Weight (kg) 17.24 kg	Colours: Black	Warranty Lifetime limited warranty on frame; 1-year limited warranty on electronics; 18-month warranty on drive control system; 2-year limited warranty on drivetrain	Notes 1. Speed varies with user weight, terrain type, battery charge, battery condition and tyre pressure. 2. Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections. 3. Battery weight may vary based on manufacturer.						