

In-Situ Recline Sling

F15

- Suitable for amputee scenarios
 - Designed to provide closed leg full-body support on transfer
 - Pommel strap for additional comfort
 - To be applied in a supine position
 - Loops supplied for provision for hip tapes
 - Soft binding as standard
- Once the patient has been transferred, this sling can be left in-situ,
 - subject to an appropriate risk assessment
 - Fixed colour-coded loop tapes as standard
 - Wash temp 71°C / 159°F
 - Tumble dry low heat (max 50°C)



NB: Silva-Safety Slot® tapes available upon request

Product Codes & Specifications:



Silva-Superfine PLUS®		A	B	C	D
XXS	F15SFPXXS	66 - 76cm	78 - 81cm	60cm	69cm
XS	F15SFPXS	76 - 86cm	81 - 92cm	70cm	75cm
S	F15SFPS	86 - 96cm	92 - 104cm	78cm	83cm
M	F15SFPM	96 - 108cm	104 - 116cm	88cm	93cm
L	F15SFPL	108 - 130cm	116 - 128cm	108cm	106cm
XL	F15SFPXL	130 - 144cm	128 - 140cm	118cm	117cm
XXL	F15SFPXXL	144 - 160cm	128 - 140cm	132cm	137cm

