

In-Situ Deluxe Leg Sling HS

S16

- Divided leg sling offering full body and head support
 - Designed to provide full body support on transfer
 - Adjustable Integral head support with removable stiffeners
 - Shaped leg sections provide full femur and hip support
 - Soft binding as standard
- Once the patient has been transferred, this sling can be left in-situ, subject to an appropriate risk assessment
 - Easy to apply with low-friction fabric to leg areas for ease of application
 - Fixed colour-coded loop tapes as standard
 - Wash temp 71°C / 159°F
 - Tumble dry low heat (max 50°C)



NB: Silva-Safety Slot® tapes available upon request

Product Codes & Specifications:

 SWL 220kg  SWL 64kg  SWL 220kg							
	Silva-Superfine PLUS®	Silva-Superfine® with COOLMAX®	Silva-Superfine® with COOLMAX® Double Layer	A	B	C	D
XXXS	S16FXXXS	S16SFCXXXSM	–	56 - 66cm	50 - 55cm	56cm	60cm
XXS	S16FXXS	S16SFCXXSM	–	66 - 76cm	56 - 62cm	64cm	65cm
XS	S16FXSM	S16SFCXSM	–	76 - 86cm	60 - 66cm	73cm	70cm
S	S16FSM	S16SFCSM	S16SFCDLS	86 - 96cm	66 - 74cm	80cm	78cm
M	S16FMED	S16SFCMED	S16SFCDLM	96 - 108cm	74 - 80cm	90cm	86cm
L	S16FLRG	–	S16SFCDLL	108 - 130cm	80 - 88cm	97cm	94cm
XL	S16FXLRG	–	S16SFCDLXL	130 - 144cm	88 - 100cm	114cm	103cm
XXL	S16FXXLRG	–	S16SFCDLXXL	144 - 160cm	96 - 104cm	120cm	107cm

